



PHOTOMULTIPLIER TUBE R12829

Times

800 nm High Sensitivity Multialkali Photocathode 28 mm (1-1/8 Inch) Diameter, 9-Stage, Side-On Type





Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult with our sales office. Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications are subject to change without notice. No patent rights are granted to any of the circuits described herein. ©2014 Hamamatsu Photonics K.K.

SPECIFICATIONS

GENERAL

	Parameter	Description / Value	Unit
Spectral Respo	onse	185 to 900	nm
Wavelength of	Maximum Response	450	nm
Photocathode	Material	Multialkali	—
Filotocathoue	Minimum Effective Area	8×24	mm
Window Mater	al	UV glass	_
Dunada	Structure	Circular-cage	—
Dynode	Number of Stages	9	—
Direct	Anode to Last Dynode	4	pF
Interelectrode	Anode to All Other	6	pF
Capacitances	Electrodes	0	μr
Base		11-pin base JEDEC No. B11-88	—
Weight		Approx. 45	g
Operating Amb	pient Temperature	-30 to +50	°C
Storage Tempe	erature	-30 to +50	°C
Suitable Socke	et	E678–11A (Sold Separately)	_
Suitable Socke	et Assembly	E717–63 (Sold Separately)	

MAXIMUM RATINGS (Absolute Maximum Values)

	Parameter	Value	Unit
Supply	Between Anode and Cathode	1250	V
Voltage	Between Anode and Last Dynode	250	V
Average Anode Current ^A		0.1	mA

CHARACTERISTICS (at 25 °C)

	Pa	ramet	er	Min.	Тур.	Max.	Unit	
	Rad	iont	at 450 nm		105		mA/W	
	nau	lanı	at 800 nm	60	65	_	mA/W	
	Qua	ntum	at 450 nm	—	29	—	%	
Cathode	Effic	iency	at 800 nm	9.3	10	—	%	
Sensitivity	Lum	inous ^B		600	650	—	μA/Im	
	Red	/ White	e Ratio ^c	—	0.45	—	—	
	Blue	e Sensit	tivity Index D	—	14	—	—	
Anode	Lum	inous ^E		1600	8500		A/Im	
Sensitivity	Lum	inious -		1000	8500	_		
Gain ^E				—	$1.3 imes10^7$	—	—	
Anode Dark Current F					2.0	10	nA	
(Supply voltage at 1×10^6 Gain)		_	2.0	10	ПА			
Time Response	Anoc	de Pulse	e Rise Time ^G	—	2.2	—	ns	
	Elec	tron Tra	ansit Time ^H	—	22 —		ns	
	Trans	it Time S	pread (T.T.S.) ¹	—	1.2	1.2 —		
Anode Cur	Anode Current Light Hysteresis		_	0.1	_	%		
Stability J Voltage Hysteresis		e Hysteresis	_	1.0	_	%		

NOTES

- A: Averaged over any interval of 30 seconds maximum.
- B: The light source is a tungsten filament lamp operated at a distribution temperature of 2856 K. Supply voltage is 100 V between the cathode and all other electrodes connected together as anode.
- C:Red / White ratio is the quotient of the cathode current measured using a red filter (Toshiba R-68) interposed between the light source and the tube by the cathode current measured with the filter removed under the same conditions as Note B.
- D: The value is cathode output current when a blue filter (Corning CS 5-58 polished to 1/2 stock thickness) is interposed between the light source and the tube under the same condition as Note B.
- E: Measured with the same light source as Note B and with the voltage distribution ratio shown in Table 1 below.


Table 1:Voltage Distribution Ratio

Electrodes	ł	<	Dy	1 C)y2	Dy3	Dy	y4	Dy	/5	Dy	6	Dy	7 [)y8	D	y9		Р
Ratio		1		1		1	1		1		1	1		1		1		1	

Supply Voltage: 1000 V, K: Cathode, Dy: Dynode, P: Anode

F: After 30 min Storage in Darkness

- G:The rise time is the time for the output pulse to rise from 10 % to 90 % of the peak amplitude when the entire photocathode is illuminated by a delta function light pulse.
- H: The electron transit time is the interval between the arrival of delta function light pulse at the entrance window of the tube and the time when the anode output reaches the peak amplitube. In measurement, the whole photocathode is illuminated.
- I: Also called transit time jitter. This is the fluctuation in electron transit time between individual pulses in the signal photoelectron mode, and may be defined as the FWHM of the frequency distribution of electron transit times
- J: Hysteresis is temporary instability in anode current after light and voltage are applied.




(1)Light Hysteresis

The tube is operated at 750 V with an anode current of 1 μ A for 5 minutes. The light is then removed from the tube for a minute. The tube is then re-illuminated by the previous light level for a minute to measure the variation.

(2)Voltage Hysteresis

The tube is operated at 300 V with an anode current of 0.1 μA for 5 minutes. The light is then removed from the tube and the supply voltage is quickly increased to 800 V. After a minute, the supply voltage is then reduced to the previous value and the tube is re-illuminated for a minute to measure the variation.





(at 1×10^6 gain, after 30 min storage in darkness)




Figure 4: Anode Luminous Sensitivity and Gain Characteristics






Figure 5: Typical Time Response


Figure 6: Typical Temperature Coefficient of Anode Sensitivity



PHOTOMULTIPLIER TUBE R12829

Figure 7: Dimensional Outline and Basing Diagram (Unit: mm)





HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater. N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Hersching am Ammersee, Germany, Telephone: (49)8152-2658 E-mail: info@hamamatsu.de France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy, Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: info@hamamatsu.de United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, United Kingdom, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk North Europe: Hamamatsu Photonics Italia S.r.L.: Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (49)02-93581741 E-mail: info@hamamatsu.it China: Hamamatsu Photonics (China) Co., Ltd.: B1201 Jiaming Center, No.27 Dongsanhuan Bellu, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com MAY 2014 IP