



®



Quantaurs-Tau


Fluorescence lifetime spectrometer C11367 series



HAMAMATSU

PHOTON IS OUR BUSINESS

Quantaurs-Tau is a compact system for measuring fluorescence lifetimes in the sub nanosecond to millisecond range. Operation is simple, just set the sample into the sample chamber, and enter a few conditions on the measurement software to measure the fluorescence lifetime and PL spectrum in a short time with high precision. In a typical measurement, analysis results are obtained in a mere 60 seconds.



Quantaurs

Tau

Fluorescence Lifetime

Fluorescence lifetime measurement

Measuring an excited-state relaxation process

The fluorescence spectrum obtained from an organic material or fluorescent probe is a vital parameter for controlling and evaluating the material functions and characteristics such as the peak wavelength and fluorescence intensity. However, a fluorescence spectrum usually shows time-integrated information, and so when the material contains multiple substances and reactive elements, their fluorescence spectrum can only be acquired as integrated information. An effective approach in such cases is to observe the light emission dynamics by making use of the time axis parameter. This is generally called fluorescence lifetime measurement, in which the time required for the substance excited by the pulsed light to return to its ground state is measured in the sub-nanosecond to millisecond region. This measurement allows obtaining more information such as multiple different fluorescence lifetimes even at the same wavelength and the percentage in which they are present within the material, etc.

Features

- High sensitivity measurement by photon counting method
- Time resolution better than 100 ps (by deconvolution)
- Cooling function for solution sample (-196 °C) (option)
- Phosphorescence measurement (option)
- Fluorescence spectrum measurement
- Space-saving, compact design

Easy and quick measurements

Emission Lifetime can be gotten easily and quickly only by putting the sample into sample box and setting the 4 measurement conditions.

7 excitation wavelength

280 nm, 340 nm, 365 nm, 405 nm, 470 nm, 590 nm, and 630 nm.

Analyzing different sample forms

Thin-film, solid, solutions and powder.

2 selection of detector

Standard


C11367-31/-34

wavelength
300 nm to 800 nm

NIR

C11367-32/-35

wavelength
380 nm to 1030 nm



Organic
metal
complex

Fluorescent
probe

Dye
sensitized
PV material


OLED
material

Quantum
dot


LED
phosphor

Fluorescence lifetime measurement is applicable to varied applications. Typical applications include electron movement and energy transfer reactions within or between organic metallic complex molecules, as well as fluorescence and phosphorescence lifetime measurement of materials essential for developing organic EL devices, FRET (fluorescence resonance energy transfer) in fluorescent proteins, and pass/fail testing of compound semiconductors for solar cells and LED, etc.


Detector
spectral
sensitivity



Software designed taking account of the measurement procedure ensures easy and quick measurements.



Multi-component fluorescence lifetime analysis and comparison



Multi-component analysis of up to five components

In fluorescence lifetime measurement, a phenomenon often occurs where the data is observed as the sum of the attenuation curves of multi-component fluorescence lifetimes. Quantaurus-Tau easily calculates the fluorescence lifetime data and component ratio of each element by using the dedicated software.


Highly accurate analysis by deconvolution

Deconvolution processing enables fluorescence lifetime analysis with high accuracy. When analyzing longer lifetime components such as phosphorescence, the "Tail Fit" function can be used instead of deconvolution processing.

Real-time display of time profiles and spectrum

Time profiles or spectrum are displayed on the monitor screen in real-time. This is a useful function for selecting the time scale during measurement or determining the analysis data range.

Multi-sample fluorescence lifetime analysis and comparison



Multiple data analysis on the same screen

Calculated fluorescence lifetime values are also displayed on the same screen for easy comparison analysis.

Comparisons under the same fitting conditions

To make comparison analysis under the same conditions, Quantaurus-Tau subjects the multiple samples to specific fitting ranges, IRF (Instrument Response Function), and parameter settings.


Graph editing with a graph setup feature

This allows you to change the range of each axis as needed on the comparison analysis screen so that the data can be edited to match your purpose. This feature also allows powerful normalizing whenever needed.

Acquired data can be easily stored as text data

The acquired data can be stored into the graph analysis software as text data by simple copy-and-paste operations.

Multi-sample PL spectrum analysis and comparison



Time-resolved spectrum display

Allows time-resolved spectrum display the greatest feature offered by streak camera systems.


Spectrum and fluorescence decay curve display

Displays the full width at half maximum (FWHM), peak position and peak intensity for each profile.



Multiple data loading and comparison on the same screen

Normalized processing makes multiple data comparison easy.


TADF of the blue OLED material



TADF (Thermally activated delayed fluorescence) is known well as the 3rd generation OLED material which is high efficiency and cost saving, furthermore can be replaced with phosphorescence material. The data shows the example of fluorescence lifetime measurement of blue TADF material. In order to achieve the high efficiency, the molecule design is the important factor to minimize the energy gap of excited singlet state (S_1) and excited triplet state (T_1).



Delayed fluorescence of TADF material results in that small gap of ΔE_{st} , which is defined as the energy gap of S_1 state and T_1 state, is in the micro second range. On the other hand, wide gap of ΔE_{st} is in milli second range.



A= PPZ-DPO, PPZ-3TPT, PPT-4TPT

D=PPZ-DPS, PXZ-DPS, DMAC-DPS


PPZ: 5-phenyl-5,10-dihydrophenazine
DPO: 2,5-diphenyl-1,3,4-oxadiazole
TPT: 3,4,5-triphenyl-1,2,4-triazole

DPS: diphenylsulphone
PXZ: phenoxazine
DMAC: 9,9-dimethyl-9,10-dihydroacridine


Data courtesy of Prof. Chihaya Adachi, Hajime Nakanotani
Center for Organic Photonics and Electronics Research, Kyushu Univ.

Q. Zhang, B. Li, S. Huang, H. Nomura, H. Tanaka and C. Adachi, *nature photonics*, 8, 326 (2014)

Temperature dependency of fluorescence lifetime with TADF material



Temperature dependence of transient PL decays of triarylboron-based OLED emitter doped in CBP (16 wt%).



The fluorescence life time measurement example to observe the temperature dependency of triarylboron-based compound which is the TADF material. Phosphorescence material as the typical OLED material has chemical behavior to decrease the ratio of the light components by the temperature increasing. As the temperature rose, TADF material, however, increased the delayed fluorescence components and resulted the delay components were activated by the heat. The temperature dependency measurement was done with the setup of Cryostat.

Data courtesy of Prof. Hironori Kaji, Atsushi Wakamiya, Katsuaki Suzuki, Institute for Chemical Research, Kyoto Univ.
Data courtesy of Prof. Chihaya Adachi, Center for Organic Photonics and Electronics Research, Kyushu Univ.

K. Suzuki, S. Kubo, K. Shizu, T. Fukushima, A. Wakamiya, Y. Murata, C. Adachi, H. Kaji, *Angew chem. Int. Ed.* 54, 15231 (2015).

Fluorescence Lifetime and Absolute PL Quantum Yield

There are two processes when substances are excited by light irradiation from the ground state to excited singlet state (S₁), then deactivated to the ground state again. One is radiative process such as fluorescence and the other is a non-radiative process released as heat.

The fluorescence lifetime τ (tau) is defined as

$$k_f + k_{nr} = 1 / \tau$$

where k_f is the radiative rate constant and k_{nr} is the non-radiative constant.

On the other hand, the PL Quantum Yield (Φ) is expressed as the ratio of the number of photons emitted from molecules (PN_{em}) to that absorbed by molecules (PN_{abs}).

$$\Phi = \text{PN}_{em} / \text{PN}_{abs}$$

The PL Quantum Yield Φ is also written as

$$\Phi = k_f / (k_f + k_{nr})$$

Thus, there is a correlation between τ (tau) and Φ as shown in the following equation, and they are very important parameters for controlling the emission mechanisms of the materials.

$$k_f = \Phi / \tau$$



Fluorescence lifetime spectrometer
C11367 series



Absolute PL quantum yield spectrometer
C11347 series





A diversified evaluation of the luminescence materials is available!

Quantaurus-Tau for measuring fluorescence lifetime and Quantaurus-QY for absolute PL quantum yield with simplified and minimized operating procedure are available for everybody.

Combination of Quantaurus-Tau and Quantaurus-QY allow users to obtain complementary analysis results.

ff luminescent characteristic of Lanthanide complex



Lanthanide compound has a characteristic of high luminescent performance and is expected to apply for the variety of functional materials such as OLED, photovoltaic or a sensor. The fluorescence lifetime of a series of Lanthanide complex was measured in the acetonitrile solution (at room temperature) and the solid state (at room temperature, at 77K). A series of Lanthanide complex has Bipyridine skeleton. The remarkable difference was observed to fluorescence lifetime depending on the temperature (T) and Quantum yield (QY) of TbIII complex which has the characteristic of thermal equilibrium caused by the energy transfer between excited triplet state of ligand and the energy level of center metal.

		Temp.	τ [ms](amp.)	QY ^a [%]
EuL	In the solid state	rt	1.27 (1.0)	52.6 (\pm 1.4)
		77 K	1.35 (1.0)	63.5 (\pm 2.7)
TbL	In acetonitrile	rt	1.55 (1.0)	12.0 (\pm 0.5)
		77 K	0.0153 (0.96) 0.00234 (0.04)	1.0 (\pm 0.2)
	In acetonitrile	77 K	1.49 (1.0)	91.5 (\pm 1.4)
		rt	n.d.	\approx 0

^a The values of Ln emission were based on the ligand excitation, and observed with Absolute PL quantum yield spectrometer C9920-02.

L: Ligand
amp.: amplitude
rt: room temperature

Data courtesy of Prof. Miki Hasegawa, Aoyama Gakuin Univ.

M. Hasegawa, H. Ohtsu, D. Kodama, T. Kasai, S. Sakurai, A. Ishii, and K. Suzuki, *New J. Chem.*, **38**, 1225 (2014)

Specifications

Type number	C11367-31	C11367-34	C11367-32	C11367-35
Sample	Solution, Thin-film	Solid (Thin-film, Powder)	Solution, Thin-film	Solid (Thin-film, Powder)
Detector type	Standard		NIR	
Wavelength range	300 nm to 800 nm		380 nm to 1030 nm	
Excitation light source	Seven types of LED light source (280 nm, 340 nm, 365 nm, 405 nm, 470 nm, 590 nm, 630 nm)			
Excitation light source switching	Software control			
Monochromator	Czerny-Turner monochromator			
Measurement time range	4 ns to 10 s / full scale			
Phosphorescence measurement	Phosphorescence excitation wavelength (280 nm, 340 nm, 365 nm, 405 nm, 442 nm, 470 nm, 589 nm, 632 nm)			
Time axis channel	512 ch, 1024 ch, 2048 ch, 4096 ch			
Total time resolution	< 1.0 ns FWHM (IRF with 365 nm LED)			
Analysis function	Fluorescence lifetime analysis (up to five components by exponential function fitting) and spectrum analysis			
Supported OS	Windows 7 (32 bit), Windows 7 (64 bit)			

Options

Sample box

● Sample box for solution sample A12178-02

A12178-02 is a sample box for measuring the solution samples (standard: compatible with 10 mm square cells) or thin film samples. The normal sample box of C11367-31 or C11367-32 is A12178-02.

● Sample box for solid sample A11551-02

A11551-02 is a sample box for measuring the powder samples or thin solid film samples. The normal sample box of C11367-34 or C11367-35 is A11551-02.

● Sample box for low temperature A11797-02

A11797-02 is a sample box for setting A11238-04 when measuring the lifetime of a solution samples at liquid nitrogen temperature.

● Sample box for cryostat Optistat DN A12268-01

Sample box for Optistat DN2 (Oxford Instruments).

Sample holder

● Sample holder for low temperature A11238-04

This is used to cool the solution sample with liquid nitrogen.

Sample case

For solution

● Side-arm cells (3 sets) A10095-02

● Sample tube for low temperature measurement (5 pcs) A10095-04

This is used to measure a sample solution at liquid nitrogen temperature.

For powder

● Laboratory dish without caps (5 pcs) A10095-01

● Laboratory dish with caps (5 sets) A10095-03

This is used for making measurements on powder samples. This contains 5 dishes made of synthetic quartz, which suppresses fluorescence and luminescence.



Light source option

● Xenon flash lamp unit for phosphorescence measurement C11567-02

Xenon flash lamp.

● Band pass filter A12991- XXX

Selectable from the wavelengths of either 280 nm or 340 nm.

● Band pass filter A13905- XXX

Selectable from among the following wavelengths: 365 nm, 405 nm, 442nm, 470 nm, 589 nm, 632 nm.

● PLP-10-XXXXTAU Laser diode head M12488 series

M12488 series are the dedicated laser diode heads for Quantaaurus-Tau, which can be used in combination with the controller and adapter.


Selectable from the following wavelength of 375 nm, 405 nm, 445 nm, 465 nm, 483 nm, 510 nm, 655 nm, 785 nm, 850 nm.

● Adapter A12487-01

A12487-01 is adapter for attaching M12488 to Quantaaurus-Tau. A12487-01 is used when excitation light source is the PLP-10.

Dimensional outlines

(unit : mm) Weight : 32 kg



Quantaaurus-QY is registered trademark of Hamamatsu Photonics K.K. (China, EU, Korea, Japan, U.S.A.)

Quantaaurus-Tau is registered trademark of Hamamatsu Photonics K.K. (China, EU, Korea, Japan, U.S.A.)

Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.

● Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.

● Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.

Specifications and external appearance are subject to change without notice.

© 2016 HAMAMATSU PHOTONICS K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.r.l.: Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-935-81-733, Fax: (39)02-935-81-741 E-mail: info@hamamatsu.it

China: Hamamatsu Photonics (China) Co., Ltd.: 1201 Tower B, Jiaming Center, 27 Dongsanhuan Beilu, Chaoyang District, 100020 Beijing, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: Hamamatsu Photonics Taiwan Co., Ltd.: 8F-3, No.158, Section2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (886)03-659-0080, Fax: (886)07-811-7238 E-mail: info@tw.hpk.co.jp

Cat. No. SHSS0011E13

APR/2016 HPK

Created in Japan