

# Multiplexing Calcium Mobilization and Membrane Potential Assays Using the FDSS6000

## Abstract

The FDSS6000 from Hamamatsu Photonics Systems uses Xenon lamps, enabling the use of fluorescent dyes ranging from UV to visible excitation wavelengths. FDSS6000 optics can read up to two excitation wavelengths and two emission wavelengths in the same assay.

Current high throughput screening (HTS) strategies use only one dye per assay, such as Fura-2-AM, Fluo-3-AM, or Fluo-4-AM (calcium mobilization), or DisBAC<sub>2</sub>(3) (membrane potential).

We present a novel multiplexing method for HTS Screening: Loading the same cells using both the UV dye Fura-2 AM and a no wash Membrane Potential Dye.

Multiplexing assays measure orthogonal cell responses (calcium mobilization and membrane potential) using one volume of valuable library compound.

## Introduction

Several groups report results using double dye loaded cells (Figure 1) for studying either calcium channels or G Protein Coupled Receptors, including Fura-2AM/DisBac<sub>2</sub>(3)<sup>1</sup>, Indo-1AM/DisBac<sub>2</sub>(3)<sup>2</sup>, Fura-2AM/DiBac<sub>4</sub>(3)<sup>3</sup>, and Fluo-4AM/Di4 ANEPPS<sup>4</sup>. These assays use single channel Spectrofluorometers for readout.

In this technical note we report the results of a multiplexed assay in a 384 well format.

### Double Parameter Fluorescent Measurement

Cells Loaded Using Both Fura-2 AM and Membrane Potential Dyes




Fig. 1: Multiplex Theory

## Materials and Methods

Adherent cells were washed and left with a 30  $\mu$ L residue. Fura-2 AM, 1.2  $\mu$ M final concentration, was incubated with the cells for 1 h, RT. Following cell washing (30  $\mu$ L residual volume) 15  $\mu$ L No Wash Membrane Potential Dye was added to 0.25 X. After 15 min incubation cells were read on the FDSS6000.

FDSS6000 optics were set to sequentially read ex 380 nm/em 540 nm for calcium mobilization (Fura-2 AM) and ex 480 nm/em 540 nm for membrane potential, in a 1.2 sec cycle time.

## Results

A recombinant cell line expressing Receptor A was used in the study. Cells were agonized against Ligand 1 and Ligand 2 (Figure 2).

Cells show calcium mobilization (decrease in fluorescence at ex 380 nm) following exposure to both Ligand 1 and 2 (Figure 2(a)). Figure 2(b) shows Ligand 1 induces membrane potential hyperpolarization while Ligand 2 induces depolarization.







Fig. 2: Fluorescence of calcium mobilization (a) and membrane potential (b). Ligands (arrow) were added at 10 sec.

Dose response curves measuring calcium mobilization and membrane hyperpolarization were built using Ligand 1 (Figure 3).



**Fig. 3:** Dose response curves expressing calcium mobilization and membrane hyperpolarization using Ligand 1 (n=6/group error bars=SEM).

Three antagonists against Receptor A were tested for potency (Figure 4). Note the difference in potency of compound Z; calcium mobilization inhibition is detected yet little membrane potential inhibition is shown. By contrast compound X inhibits both calcium mobilization and membrane potential with equivalent potency.



**Fig. 4:** Rank order potency of three compound antagonists using Ligand 1 at 168  $\mu\text{M}$ , (n=4/group, error bars=SEM).

## Summary

Multiplexing, simultaneously measuring changes in calcium mobilization and membrane potential, is the next step in the evolution of cell based HTS. This report shows how multiplexing further categorizes compounds, i.e. concomitant inhibition of calcium mobilization and membrane potential versus calcium mobilization inhibition only.

## References

- Montero, M. J. Garcia-Sancho, J. Alvarez. Activation by Chemotactic Peptide of a Receptor-operated  $\text{Ca}^{2+}$  entry Pathway in differentiated HL60 cells. *J. Biol. Chem.* 269: 29451-29456. 1994.
- Kremer, S.G., W. Zeng, K. Skorecki. Simultaneous fluorescence measurement of calcium and membrane potential responses to endothelin. *Am. J. Physiol.* 263 (Cell Physiol. 32): C1302-C1309. 1992.
- Laskey, R.E., D.J. Adams, M. Cannell, C. van Breeman. Calcium entry-dependent oscillations of cytoplasmic calcium concentration in cultured endothelial cell monolayers. *Proc. Natl. Acad. Sci.* 89:1690-1694. 1992.
- Johnson, P.L., W. Smith, T.C. Baynham, S.B. Knisley. Errors caused by combination of Di-4 ANEPs and Fluo 3 / 4 for simultaneous measurements of transmembrane potentials and intracellular calcium. *Ann. Biomed. Eng.* 27:563-571. 1999.

★ FDSS is registered trademark of Hamamatsu Photonics K.K. (China, France, Germany, Italy, Japan, U.K., U.S.A.)

★ Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.

- Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearance are subject to change without notice.

© 2013 Hamamatsu Photonics K.K.

**HAMAMATSU PHOTONICS K.K.** [www.hamamatsu.com](http://www.hamamatsu.com)

**HAMAMATSU PHOTONICS K.K., Systems Division**

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: [export@sys.hpk.co.jp](mailto:export@sys.hpk.co.jp)

U.S.A.: Hamamatsu Corporation, 360 Foothill Road, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: [usa@hamamatsu.com](mailto:usa@hamamatsu.com)

Germany: Hamamatsu Photonics Deutschland GmbH, Arzbergerstr. 10, D-82211 Hensching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658 E-mail: [info@hamamatsu.de](mailto:info@hamamatsu.de)

France: Hamamatsu Photonics France S.A.R.L., 19, Rue du Saule Trépas, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: [infos@hamamatsu.fr](mailto:infos@hamamatsu.fr)

United Kingdom: Hamamatsu Photonics UK Limited, 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, United Kingdom, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: [info@hamamatsu.co.uk](mailto:info@hamamatsu.co.uk)

North Europe: Hamamatsu Photonics Norden AB, Torshamnsgatan 35 SE-164 40 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: [info@hamamatsu.se](mailto:info@hamamatsu.se)

Italy: Hamamatsu Photonics Italia S.r.l., Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-93581733, Fax: (39)02-93581741 E-mail: [info@hamamatsu.it](mailto:info@hamamatsu.it)

China: Hamamatsu Photonics (China) Co., Ltd., B1201 Jiaming Center, No.27 Dongshanhuan Bellu, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: [hpc@hamamatsu.com.cn](mailto:hpc@hamamatsu.com.cn)

Cat. No. SBI0028E02

NOV/2013 IP

Created in Japan (PDF)