

# EFS pacing kit for FDSS/μCELL


Specialized for Cardiomyocyte Assays

## Overview

- Synchronizing cardiomyocyte beating in all 96 wells in a plate by electric field stimulation (EFS)
- Pacing cardiomyocytes with variable frequencies, with variable voltage and pulse duration of electric field stimulation

Electric Field Stimulation (EFS) pacing kit for FDSS/μCELL is an option to have an ability to stimulate cardiomyocytes electrically by the 96-electrode array to pace beatings of the cells, in which the cells are cultivated in a 96-well plate format. This option could be associated with other options available in FDSS/μCELL, such as assay plate heating unit and fast data acquisition software (up to 120 Hz frame reading), each of which could make more accurate monitoring of intracellular molecular events such as  $\text{Ca}^{2+}$  transients in cardiomyocytes. The EFS option with FDSS/μCELL should be useful in the drug discovery as in cardiomyocyte research.

## Rat primary cardiomyocytes in 96-well plate, intracellular $\text{Ca}^{2+}$ transient (Fluo-8)



### Caution Notice:


The FDSS/μCELL EFS system should not be used for optically detecting/monitoring change in transmembrane potential of the cells.

The FDSS/μCELL EFS system should not be used on any cell or cells in which the user or anyone else has expressed target ion channels.

# FDSS Option

## Measurement Example

### Ca<sup>2+</sup> transient measurement with EFS stimulation at various frequencies



▲ FDSS/μCELL

- Fluorescence / Luminescence measurements while electric field stimulation.
- Variable EFS frequency.
- Variable EFS voltage.

\* EFS pacing kit option requires fluorescence / luminescence sensor and the latest data analysis computer.  
Please contact your local sales representative for more detail.

Part number	Description
A11529-14	Heater Option
U8524-11	High Speed Acquisition Software
U8524-12	Analysis software for cardiomyocyte
M13040-01	96ch EFS Pacing Stimulation kit

- ★ FDSS is registered trademark of Hamamatsu Photonics K.K. (China, France, Germany, Italy, Japan, U.K., U.S.A.)
- ★ Hamamatsu Photonics K.K. and its subsidiaries or representatives do not warrant that the use of our products and those goods does not infringe the intellectual property rights of third parties.
- ★ Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.
- Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearance are subject to change without notice.

© 2014 Hamamatsu Photonics K.K.

**HAMAMATSU PHOTONICS K.K.** [www.hamamatsu.com](http://www.hamamatsu.com)

**HAMAMATSU PHOTONICS K.K., Systems Division**

**812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp**

**U.S.A.:** Hamamatsu Corporation; 360 Foothill Road, Bridgewater, N.J 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

**Germany:** Hamamatsu Photonics Deutschland GmbH; Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

**France:** Hamamatsu Photonics France S.A.R.L.; 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

**United Kingdom:** Hamamatsu Photonics UK Limited; 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

**North Europe:** Hamamatsu Photonics Norden AB; Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

**Italy:** Hamamatsu Photonics Italia S.r.l.; Strada della Moia, 1 int. 6 20020 Arese (Milano), Italy, Telephone: (39)02-93581733, Fax: (39)02-93581741 E-mail: info@hamamatsu.it

**China:** Hamamatsu Photonics (China) Co., Ltd.; B1201 Jiaming Center, No.27 Dongsanhuan Beilu, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Cat. No. SBIS0106E02  
NOV/2014 HPK  
Created in Japan