

Photon counting module



C14463-050GD

Fiber coupling type, low-light-level detection

The C14463-050GD is a photon counting module that can detect low-level light. It consists of a TE-cooled single pixel photon counter (SPPC), amplifier, comparator, high-voltage power supply circuit, and temperature control circuit. This product is a fiber coupling type module. You can simply supply external power (± 5 V) to use this module.

Features

- ➡ Fiber coupling type
- **■** Built-in single pixel photon counter
- High sensitivity in the visible to near IR region
- **■** Low dark count

Applications

- **Low-light-level measurement**
- **▶** Particle diameter measurement
- **➡** Fluorescence measurement
- Analytical instruments

- Absolute maximum ratings

Parameter	Symbol	Condition	Value	Unit
Supply voltage	Vs		±6	V
Operating temperature	Topr	No dew condensation*1	-10 to +40	°C
Storage temperature	Tstg	No dew condensation*1	-20 to +70	°C

^{*1:} When there is a temperature difference between a product and the surrounding area in high humidity environment, dew condensation may occur on the product surface. Dew condensation on the product may cause deterioration in characteristics and reliability.

- Recommended operating conditions

Parameter		Symbol	Min.	Тур.	Max.	Unit
Supply voltage*2	Positive power supply	Vs	4.75	5	5.25	V
	Negative power supply		-4.75	-5	-5.25	

^{*2:} A power supply with 1 A or higher output must be used.


Electrical and optical characteristics (Ta=25 °C, λ=600 nm, Vs=±5 V, unless otherwise noted)

Paramete	r	Symbol	Condition	Min.	Typ.	Max.	Unit
Spectral response range		λ		370 to 1000			nm
Peak sensitivity waveleng	gth	λр		- 600 -		-	nm
Fiber connector*3 -			FC type			-	
Chip temperature (settin	g temperature)	Tchip		20 -		°C	
Photon detection efficiency F		PDE		25	35	-	%
Dark count		CD		-	20	60	cps
Afterpulse probability		-	100 ns to 500 ns	- 0.1 -		%	
Comparator output	ator output - TTL compatible			-			
Current consumption	+5 V	Ic		-	+200	+1000	mA
	-5 V			-	-20	-40	

^{*3:} Recommended fiber: GI 50/125 multimode fiber


Note: Exceeding the absolute maximum ratings even momentarily may cause a drop in product quality. Always be sure to use the product within the absolute maximum ratings.

Block diagram




Connection example


Using the supplied power cable, connect the module to a power supply. You can count output pulses by connecting the module to a frequency counter.




Measurement example




- Photon detection efficiency vs.wavelength (typical example)



Linearity (typical example)



- Dimensional outline (unit: mm)



KACCA0394EA

Accessories

- · Power cable
- · Instruction manual

Options (sold separately)

Coaxial conversion adapter A10613 series

These are coaxial conversion adapters for converting the SMB coaxial connector for extracting MPPC module signals into a BNC coaxial connector or an SMA coaxial connector. These adapters make connection to a BNC cable or SMA cable possible.





A10613-01 (SMB-BNC)

A10613-02 (SMB-SMA)

Precautions

· Use the product by referring to the supplied instruction manual.

- Related information

http://www.hamamatsu.com/sp/ssd/doc_en.html

- Precautions
- Disclaimer

Information described in this material is current as of March 2020.

Product specifications are subject to change without prior notice due to improvements or other reasons. This document has been carefully prepared and the information contained is believed to be accurate. In rare cases, however, there may be inaccuracies such as text errors. Before using these products, always contact us for the delivery specification sheet to check the latest specifications.

The product warranty is valid for one year after delivery and is limited to product repair or replacement for defects discovered and reported to us within that one year period. However, even if within the warranty period we accept absolutely no liability for any loss caused by natural disasters or improper product use. Copying or reprinting the contents described in this material in whole or in part is prohibited without our prior permission.

MAMATSU

www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Solid State Division

1126-1 Ichino-cho, Higashi-ku, Hamamatsu City, 435-8558 Japan, Telephone: (81)53-434-3311, Fax: (81)53-434-5184

U.S.A.: Hamamatsu Corporation: 360 Footbill Road, Bridgewater, N.J. 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218, E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8, E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy, Cedex, France, Telephone: (43)16 95 37 1 00, Fax: (33)1 69 53 71 10, E-mail: info@hamamatsu.fr

United Kingdom: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01, E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.r.l.: Strada della Moia, 1 int. 6, 20020 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41, E-mail: info@hamamatsu.it

China: Hamamatsu Photonics (China) Co., Ltd.: B1201, Jiaming Center, No.27 Dongsanhuan Bellu, Chaoyang District, 100020 Beijing, P.R.China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866, E-mail: hpc@hamamatsu.com.cn

Taiwan: Hamamatsu Photonics Taiwan Co., Ltd.: 8F-3, No. 158, Section2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C. Telephone: (86)3-659-0081, E-mail: info@hamamatsu.com.cn